

6000  kg



BIMI

SMART BODY SCALE

User Manual

PRODUCT INFORMATION

Biometric Measurement

Weight, BMI, BMR.

Measurement Accuracy

50KG: ±100g	100KG: ±200g	150KG: ±300g
100lbs: ±0.22lbs	220lbs: ±0.44lbs	331lbs: ±0.66lbs

PRODUCT SPECIFICATIONS

- **Device available in:** white, black
- **Weight Unit:** kg/lb
- **Dimensions:** 280*280*25mm/11.02*11.02*0.98in
- **LED Screen:** 63.5*26mm / 2.5*1.02in
- **Maximum Weight Range:** 180kg/396lb
- **Minimum Weight Range:** 0.1kg/0.2lb
- **Operating weight:** ≥5kg
- **Working Temperature Range:** 5~40°C
- **Working Humidity Range:** ≤90%RH
- **Storage Temperature Range:** -10°C~60°C
- **Batteries:** 3 AAA batteries (included)

QUICK START WITH YOUR SMARTPHONE


Step 1: Download the NOERDEN App

To download the App, scan the QR code or search **NOERDEN** in the App Store or on Google Play and follow the instructions.




Step 2: Connect Your Device

Insert the supplied batteries in your device and turn your phone's Bluetooth ON. Then go to the "**Device**" page in the App to add your new device and follow the instructions.



Step 3: Analyze Your Body Composition

Step on the body scale with bare feet. Then go to the “**Weigh-in**” section of the homepage on the App to start the analysis. To access your detailed data just click on this page.



Step 4: How to check previous body analysis measurements

Click on the “**Weigh-in**” section of the **NOERDEN** App homepage. You will be redirected to your Body Composition detailed data.



CONDITIONS OF USE

- **Storage & Usage Conditions:** use the scale on a flat, hard surface. Avoid direct sunlight and humid environments. Keep at temperature range -10°C to 60°C (14°F to 140°F.)
- **Minimum recommended age for product use is 10 years old.** However, the minimum age by law to create a user account is 15 years old. In any case, monitor children while they use **NOERDEN** products.
- **Battery:** to change the batteries, open the compartment located on the bottom of the body scale. Please change the batteries when you see the "LO" indicator displayed on the scale's screen.
- **Automatic Shutdown:** the device will automatically shutdown if (1) 10 seconds have passed after its last use, (2) if the maximum load is exceeded by 2.5kg/5.5lb/0.4st

Compatibility

Please check your phone compatibility before using the device:

Compatible with iPhone 6 and iOS 12.0 or above;

Compatible with smartphones with Android 6.0 or above;

Bluetooth 4.0 .

FAQ

Q: Why doesn't my scale work?/Why doesn't my scale display anything?

A: Please make sure the batteries are properly inserted, the film covering the display has been removed, and the scale is not in low-battery mode.

Q: Why does the screen say "LO" when I turn on my scale?

A: The batteries are low, please replace them.

Q: Why do the results differ based on time and place?

A: Body composition changes constantly and may slightly differ during different times. Make sure that you take measurements under the same conditions for accurate results. Also, please remember to place the scale on a hard, flat surface.

Q: Why does the screen display "ERR"?

A: Reset the scale by removing and replacing the batteries. Please don't step on the scale until the screen displays "0.0".

ATTENTION

- These results are for reference only. Do not use it as a basis for your diet and fitness regiment.
- Do not hop or jump onto the scale. Always stand with both feet flat in the center.
- To avoid injury, do not allow children to play with the scale.
- Please insert the batteries in the correct position and change them regularly.
- To avoid damage, please make sure your feet are dry and clean before use.
- For the most accurate results, please take your measurement at the same time every day.
- Do not use the scale after strenuous exercise.
- Do not disassemble the scale by yourself.
- Do not let the scale come into contact with water.
- Do not use hot water and harsh detergents to clean the scale.

DECLARATION OF CONFORMITY

The undersigned, **NOERDEN** declares that the electrical radio equipment types of the **NOERDEN** Smart Body Scales, Model BIM1, comply with the 2014/53/EU directive.

www.noerden.io

For more information,
go to the FAQ page in the App or on our website